



X Antibody, Biotin conjugated

Product Code	CSB-PA323774LD01WAW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17401
Immunogen	Recombinant Woodchuck hepatitis B virus Protein X protein (1-141AA)
Raised In	Rabbit
Species Reactivity	Woodchuck hepatitis B virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein X (HBx) (Peptide X) (pX), X
Immunogen Species	Woodchuck hepatitis B virus (isolate 8) (WHV)
Research Area	Others
Target Names	X
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.