



# vpx Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA323941HC01SHD
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P19508
<b>Immunogen</b>	Recombinant Simian immunodeficiency virus Protein Vpx protein (1-112AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Simian immunodeficiency virus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein Vpx (Viral protein X) (X ORF protein), vpx
<b>Immunogen Species</b>	Simian immunodeficiency virus (isolate PBj14/BCL-3) (SIV-sm) (Simian immunodeficiency virus sooty mangabey monkey)
<b>Research Area</b>	Others
<b>Target Names</b>	vpx
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.