



vpx Antibody, Biotin conjugated

Product Code	CSB-PA323941HD01SHD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19508
Immunogen	Recombinant Simian immunodeficiency virus Protein Vpx protein (1-112AA)
Raised In	Rabbit
Species Reactivity	Simian immunodeficiency virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Vpx (Viral protein X) (X ORF protein), vpx
Immunogen Species	Simian immunodeficiency virus (isolate PBJ14/BCL-3) (SIV-sm) (Simian immunodeficiency virus sooty mangabey monkey)
Research Area	Others
Target Names	vpx
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.