



toxA Antibody, Biotin conjugated

Product Code	CSB-PA324524HD01DUN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16154
Immunogen	Recombinant Peptoclostridium difficile Toxin A protein (2387-2710AA)
Raised In	Rabbit
Species Reactivity	Peptoclostridium difficile
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin A (EC 3.4.22.-), toxA, tcdA
Immunogen Species	Clostridioides difficile (Peptoclostridium difficile)
Research Area	Others
Target Names	toxA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.