



P Antibody, HRP conjugated

Product Code	CSB-PA324598LB01HDO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17100
Immunogen	Recombinant Hepatitis B virus genotype A2 subtype adw2 Protein P protein (349-692AA)
Raised In	Rabbit
Species Reactivity	Hepatitis B virus genotype A2 subtype adw2
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein P [Includes: DNA-directed DNA polymerase (EC 2.7.7.7); RNA-directed DNA polymerase (EC 2.7.7.49); Ribonuclease H (EC 3.1.26.4)], P
Immunogen Species	Hepatitis B virus genotype A2 subtype adw2 (isolate Germany/991/1990) (HBV-A)
Research Area	Others
Target Names	P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.