



# P Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA324598LD01HDO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17100
<b>Immunogen</b>	Recombinant Hepatitis B virus genotype A2 subtype adw2 Protein P protein (349-692AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Hepatitis B virus genotype A2 subtype adw2
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein P [Includes: DNA-directed DNA polymerase (EC 2.7.7.7); RNA-directed DNA polymerase (EC 2.7.7.49); Ribonuclease H (EC 3.1.26.4)], P
<b>Immunogen Species</b>	Hepatitis B virus genotype A2 subtype adw2 (isolate Germany/991/1990) (HBV-A)
<b>Research Area</b>	Others
<b>Target Names</b>	P
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.