



KP6 killer toxin Antibody, HRP conjugated

Product Code	CSB-PA325875LB01UBC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16948
Immunogen	Recombinant Ustilago maydis P6 virus KP6 killer toxin protein (28-105AA)
Raised In	Rabbit
Species Reactivity	Ustilago maydis P6 virus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KP6 killer toxin (Killer protein 6) [Cleaved into: KP6 killer toxin subunit alpha (VP10); KP6 killer toxin subunit beta (VP12.5)]
Immunogen Species	Ustilago maydis P6 virus (UmV6) (UmV-P6)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.