



Ribosome-inactivating protein karasurin-C Antibody, Biotin conjugated

Product Code	CSB-PA326146HD01TIF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24478
Immunogen	Recombinant Trichosanthes kirilowii Ribosome-inactivating protein karasurin-C protein (22-270AA)
Raised In	Rabbit
Species Reactivity	Trichosanthes kirilowii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein karasurin-C (EC 3.2.2.22) (rRNA N-glycosidase) [Cleaved into: Ribosome-inactivating protein karasurin-A]
Immunogen Species	Trichosanthes kirilowii (Chinese snake gourd) (Chinese cucumber)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.