



EPB2 Antibody, Biotin conjugated

Product Code	CSB-PA326194LD01HWQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25250
Immunogen	Recombinant Hordeum vulgare Cysteine proteinase EP-B 2 protein (134-373AA)
Raised In	Rabbit
Species Reactivity	Hordeum vulgare
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cysteine proteinase EP-B 2 (EC 3.4.22.-), EPB2
Immunogen Species	Hordeum vulgare (Barley)
Research Area	Others
Target Names	EPB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.