



# AAP1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA327747LC01SVG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P37898
<b>Immunogen</b>	Recombinant Saccharomyces cerevisiae Alanine/arginine aminopeptidase protein (1-389AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Saccharomyces cerevisiae
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Alanine/arginine aminopeptidase (EC 3.4.11.-), AAP1
<b>Immunogen Species</b>	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Research Area</b>	Others
<b>Target Names</b>	AAP1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.