



Nuclease P1 Antibody, Biotin conjugated

Product Code	CSB-PA328459LD01PEU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24289
Immunogen	Recombinant Penicillium citrinum Nuclease P1 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Penicillium citrinum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclease P1 (EC 3.1.30.1) (Deoxyribonuclease P1) (Endonuclease P1)
Immunogen Species	Penicillium citrinum
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.