



cfaE Antibody

Product Code	CSB-PA328563LA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25734
Immunogen	Recombinant Escherichia coli CFA/I fimbrial subunit E protein (1-360AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CFA/I fimbrial subunit E (Colonization factor antigen I subunit E), cfaE
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	cfaE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.