



tcyA Antibody, Biotin conjugated

Product Code	CSB-PA331446HD01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42199
Immunogen	Recombinant Bacillus subtilis L-cystine-binding protein TcyA protein (20-268AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-cystine-binding protein TcyA, tcyA, yckK
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	tcyA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.