



# GV51\_0914 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA333053LC01GHL   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                       |
| <b>Uniprot No.</b>         | D6SZ78  |
| <b>Immunogen</b>           | Recombinant Gardnerella vaginalis 5-1 Modification methylase Accl protein (1-525AA) |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Gardnerella vaginalis 5-1   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4    |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | GV51_0914   |
| <b>Immunogen Species</b>   | Gardnerella vaginalis 5-1   |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | GV51_0914   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.         |