



GV51_0914 Antibody, Biotin conjugated

Product Code	CSB-PA333053LD01GHL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D6SZ78
Immunogen	Recombinant Gardnerella vaginalis 5-1 Modification methylase Accl protein (1-525AA)
Raised In	Rabbit
Species Reactivity	Gardnerella vaginalis 5-1
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GV51_0914
Immunogen Species	Gardnerella vaginalis 5-1
Research Area	Others
Target Names	GV51_0914
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.