



aac Antibody, HRP conjugated

Product Code	CSB-PA333489LB01ACY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29958
Immunogen	Recombinant Actinoplanes utahensis Aculeacin-A acylase protein (35-214AA)
Raised In	Rabbit
Species Reactivity	Actinoplanes utahensis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aculeacin-A acylase (EC 3.5.1.-) [Cleaved into: Aculeacin-A acylase small subunit; Aculeacin-A acylase large subunit], aac
Immunogen Species	Actinoplanes utahensis
Research Area	Others
Target Names	aac
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.