



R107.2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA333729LA01CXY |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P32740 |
| Immunogen | Recombinant Caenorhabditis elegans ATP-dependent (S)-NAD(P)H-hydrate dehydratase protein (1-307AA) |
| Raised In | Rabbit |
| Species Reactivity | Caenorhabditis elegans |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | ATP-dependent (S)-NAD(P)H-hydrate dehydratase (EC 4.2.1.93) (ATP-dependent NAD(P)HX dehydratase) |
| Immunogen Species | Caenorhabditis elegans |
| Research Area | Others |
| Target Names | R107.2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |