



phyA Antibody

Product Code	CSB-PA333934LA01DOZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34752
Immunogen	Recombinant Aspergillus niger 3-phytase A protein (24-467AA)
Raised In	Rabbit
Species Reactivity	Aspergillus niger
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-phytase A (EC 3.1.3.8) (3 phytase A) (Myo-inositol hexakisphosphate phosphohydrolase A) (Myo-inositol-hexaphosphate 3-phosphohydrolase A), phyA
Immunogen Species	Aspergillus niger
Research Area	Others
Target Names	phyA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.