



phyA Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA333934LB01DOZ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P34752 |
| Immunogen | Recombinant Aspergillus niger 3-phytase A protein (24-467AA) |
| Raised In | Rabbit |
| Species Reactivity | Aspergillus niger |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 3-phytase A (EC 3.1.3.8) (3 phytase A) (Myo-inositol hexakisphosphate phosphohydrolase A) (Myo-inositol-hexaphosphate 3-phosphohydrolase A), phyA |
| Immunogen Species | Aspergillus niger |
| Research Area | Others |
| Target Names | phyA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |