



E6 Antibody, HRP conjugated

Product Code	CSB-PA334186LB01HNY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36814
Immunogen	Recombinant Human papillomavirus type 52 Protein E6 protein (1-148AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 52
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E6, E6
Immunogen Species	Human papillomavirus type 52
Research Area	Others
Target Names	E6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.