



SUS1 Antibody, HRP conjugated

Product Code	CSB-PA335680LB01SVV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31925
Immunogen	Recombinant Saccharum officinarum Sucrose synthase protein (1-218AA)
Raised In	Rabbit
Species Reactivity	Saccharum officinarum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sucrose synthase (EC 2.4.1.13) (Sucrose-UDP glucosyltransferase) (Fragment), SUS1
Immunogen Species	Saccharum officinarum (Sugarcane)
Research Area	Others
Target Names	SUS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.