



PPDK Antibody, HRP conjugated

Product Code	CSB-PA336270LB01EKM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37213
Immunogen	Recombinant Entamoeba histolytica Pyruvate, phosphate dikinase protein (1-342AA)
Raised In	Rabbit
Species Reactivity	Entamoeba histolytica
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate, phosphate dikinase (EC 2.7.9.1) (Pyruvate, orthophosphate dikinase), PPDK
Immunogen Species	Entamoeba histolytica
Research Area	Others
Target Names	PPDK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.