



PPDK Antibody, Biotin conjugated

Product Code	CSB-PA336270LD01EKM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37213
Immunogen	Recombinant Entamoeba histolytica Pyruvate, phosphate dikinase protein (1-342AA)
Raised In	Rabbit
Species Reactivity	Entamoeba histolytica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate, phosphate dikinase (EC 2.7.9.1) (Pyruvate, orthophosphate dikinase), PPDK
Immunogen Species	Entamoeba histolytica
Research Area	Others
Target Names	PPDK