



# PPDK Antibody, Biotin conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA336270LD01EKM   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                       |
| <b>Uniprot No.</b>         | P37213  |
| <b>Immunogen</b>           | Recombinant Entamoeba histolytica Pyruvate, phosphate dikinase protein (1-342AA)    |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Entamoeba histolytica   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Biotin  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4    |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Pyruvate, phosphate dikinase (EC 2.7.9.1) (Pyruvate, orthophosphate dikinase), PPDK |
| <b>Immunogen Species</b>   | Entamoeba histolytica   |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | PPDK  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.         |