



CYCB1-1 Antibody, FITC conjugated

Product Code	CSB-PA339075LC01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30183
Immunogen	Recombinant Arabidopsis thaliana Cyclin-B1-1 protein (1-428AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-B1-1 (Cyc1-At) (G2/mitotic-specific cyclin-B1-1) (CycB1;1), CYCB1-1, CYC1
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	CYCB1-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.