



Ribosome-inactivating protein bryodin I Antibody, HRP conjugated

Product Code	CSB-PA339327LB01BMR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33185
Immunogen	Recombinant Bryonia dioica Ribosome-inactivating protein bryodin I protein (24-247AA)
Raised In	Rabbit
Species Reactivity	Bryonia dioica
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein bryodin I (EC 3.2.2.22) (BD1) (rRNA N-glycosidase)
Immunogen Species	Bryonia dioica (Red bryony) (Bryonia cretica subsp. dioica)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.