



entC2 Antibody, HRP conjugated

Product Code	CSB-PA339410HB01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34071
Immunogen	Recombinant Staphylococcus aureus Enterotoxin type C-2 protein (28-266AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enterotoxin type C-2 (SEC2), entC2
Immunogen Species	Staphylococcus aureus
Research Area	Microbiology
Target Names	entC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.