



cwIM Antibody

Product Code	CSB-PA339780HA01BQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37134
Immunogen	Recombinant Bacillus licheniformis N-acetylmuramoyl-L-alanine amidase CwIM protein (1-253AA)
Raised In	Rabbit
Species Reactivity	Bacillus licheniformis
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylmuramoyl-L-alanine amidase CwIM (EC 3.5.1.28) (Autolysin) (Cell wall hydrolase), cwIM
Immunogen Species	Bacillus licheniformis
Research Area	Others
Target Names	cwIM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.