



cwIM Antibody, HRP conjugated

Product Code	CSB-PA339780HB01BQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37134
Immunogen	Recombinant Bacillus licheniformis N-acetylmuramoyl-L-alanine amidase CwIM protein (1-253AA)
Raised In	Rabbit
Species Reactivity	Bacillus licheniformis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylmuramoyl-L-alanine amidase CwIM (EC 3.5.1.28) (Autolysin) (Cell wall hydrolase), cwIM
Immunogen Species	Bacillus licheniformis
Research Area	Others
Target Names	cwIM