



pal Antibody, HRP conjugated

Product Code	CSB-PA340792HB01LNY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26493
Immunogen	Recombinant Legionella pneumophila Peptidoglycan-associated lipoprotein protein (22-176AA)
Raised In	Rabbit
Species Reactivity	Legionella pneumophila
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan-associated lipoprotein (19 kDa surface antigen) (PPL), pal, ppIA
Immunogen Species	Legionella pneumophila
Research Area	Others
Target Names	pal
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.