



pal Antibody, HRP conjugated

Product Code	CSB-PA340792HB01LNY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26493
Immunogen	Recombinant Legionella pneumophila Peptidoglycan-associated lipoprotein protein (22-176AA)
Raised In	Rabbit
Species Reactivity	Legionella pneumophila
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan-associated lipoprotein (19 kDa surface antigen) (PPL), pal, pplA
Immunogen Species	Legionella pneumophila
Research Area	Others
Target Names	pal