




IFNAG Antibody

Product Code	CSB-PA344226LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49877
Immunogen	Recombinant Bovine Interferon alpha-G protein (24-189AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-G (IFN-alpha7), IFNAG, IFNA7
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	IFNAG

Image



Western blot
 All lanes: IFNAG antibody at 2µg/ml + HL60 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 22 kDa
 Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.