



vacA Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA345927LC01B2M |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P55981 |
| Immunogen | Recombinant Helicobacter pylori Vacuolating cytotoxin autotransporter protein (37-245AA) |
| Raised In | Rabbit |
| Species Reactivity | Helicobacter pylori |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Vacuolating cytotoxin autotransporter [Cleaved into: Vacuolating cytotoxin; Vacuolating cytotoxin translocator], vacA |
| Immunogen Species | Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori) |
| Research Area | Others |
| Target Names | vacA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |