



vacA Antibody, Biotin conjugated

Product Code	CSB-PA345927LD01B2M
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55981
Immunogen	Recombinant Helicobacter pylori Vacuolating cytotoxin autotransporter protein (37-245AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolating cytotoxin autotransporter [Cleaved into: Vacuolating cytotoxin; Vacuolating cytotoxin translocator], vacA
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Research Area	Others
Target Names	vacA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.