



POLR2J Antibody, Biotin conjugated

Product Code	CSB-PA347112LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52435
Immunogen	Recombinant Human DNA-directed RNA polymerase II subunit RPB11-a protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase II subunit RPB11-a (RNA polymerase II subunit B11-a) (RPB11a) (DNA-directed RNA polymerase II subunit J-1) (RNA polymerase II 13.3 kDa subunit), POLR2J, POLR2J1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLR2J
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.