



Delta-ctenitoxin-Pn1a Antibody

Product Code	CSB-PA349440HA01EUV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59368
Immunogen	Recombinant Phoneutria nigriventer Delta-ctenitoxin-Pn1a protein (35-82AA)
Raised In	Rabbit
Species Reactivity	Phoneutria nigriventer
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delta-ctenitoxin-Pn1a (Delta-CNTX-Pn1a) (Insecticidal neurotoxin Tx4(6-1)) (PNTx4(6-1))
Immunogen Species	Phoneutria nigriventer (Brazilian armed spider) (Ctenus nigriventer)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.