




CMC4 Antibody

Product Code	CSB-PA349766LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56277
Immunogen	Recombinant Human Cx9C motif-containing protein 4 protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cx9C motif-containing protein 4 (Mature T-cell proliferation 1 neighbor protein) (Mature T-cell proliferation-1 type A) (MTCP-1 type A) (Protein p8 MTCP-1) (p8MTCP1), CMC4, C6.1B MTCP1 MTCP1NB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CMC4

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA349766LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.