



fkpA Antibody, HRP conjugated

Product Code	CSB-PA355178HB01EGX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P65764
Immunogen	Recombinant Escherichia coli O6:H1 FKBP-type peptidyl-prolyl cis-trans isomerase FkpA protein (26-270AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O6:H1
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FKBP-type peptidyl-prolyl cis-trans isomerase FkpA (PPIase) (EC 5.2.1.8) (Rotamase), fkpA
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Research Area	Others
Target Names	fkpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.