



# gag Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA356017EC01MHK
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P03332
<b>Immunogen</b>	Recombinant Moloney murine leukemia virus Gag polyprotein protein (216-478AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Moloney murine leukemia virus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Gag polyprotein (Pr65gag) (Core polyprotein) [Cleaved into: Matrix protein p15 (MA); RNA-binding phosphoprotein p12 (pp12); Capsid protein p30 (CA); Nucleocapsid protein p10-Gag (NC-gag)], gag
<b>Immunogen Species</b>	Moloney murine leukemia virus
<b>Research Area</b>	Others
<b>Target Names</b>	gag
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.