



# F Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA356041LB01HTO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P03420
<b>Immunogen</b>	Recombinant Human respiratory syncytial virus A Fusion glycoprotein F0 protein (27-529AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human respiRatory syncytial virus A
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fusion glycoprotein F0 (Protein F) [Cleaved into: Fusion glycoprotein F2'; Interchain peptide; Fusion glycoprotein F2; Fusion glycoprotein F1], F
<b>Immunogen Species</b>	Human respiratory syncytial virus A (strain A2)
<b>Research Area</b>	Others
<b>Target Names</b>	F
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.