



F Antibody, Biotin conjugated

Product Code	CSB-PA356041LD01HTO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03420
Immunogen	Recombinant Human respiratory syncytial virus A Fusion glycoprotein F0 protein (27-529AA)
Raised In	Rabbit
Species Reactivity	Human respiRatory syncytial virus A
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fusion glycoprotein F0 (Protein F) [Cleaved into: Fusion glycoprotein F2'; Interchain peptide; Fusion glycoprotein F2; Fusion glycoprotein F1], F
Immunogen Species	Human respiratory syncytial virus A (strain A2)
Research Area	Others
Target Names	F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.