



UL48 Antibody, Biotin conjugated

Product Code	CSB-PA356268HD01HWZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04486
Immunogen	Recombinant Human herpesvirus 1 Tegument protein VP16 protein (1-490AA)
Raised In	Rabbit
Species Reactivity	Human herpesvirus 1
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tegument protein VP16 (Alpha trans-inducing protein) (Alpha-TIF) (ICP25) (Vmw65)
Immunogen Species	Human herpesvirus 1 (strain F) (HHV-1) (Human herpes simplex virus 1)
Research Area	Others
Target Names	UL48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.