



HA Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA356317LB01IDG |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04664 |
| Immunogen | Recombinant Influenza A virus Hemagglutinin protein (1-328AA) |
| Raised In | Rabbit |
| Species Reactivity | Influenza A virus |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain] (Fragment), HA |
| Immunogen Species | Influenza A virus (strain A/England/878/1969 H3N2) |
| Research Area | Others |
| Target Names | HA |