



Histidine-rich protein PFHRP-II Antibody, HRP conjugated

Product Code	CSB-PA356476HB01PLT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05227
Immunogen	Recombinant Plasmodium falciparum Histidine-rich protein PFHRP-II protein (31-289AA)
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histidine-rich protein PFHRP-II
Immunogen Species	Plasmodium falciparum
Research Area	Others