



ompC Antibody, HRP conjugated

Product Code	CSB-PA356994LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06996
Immunogen	Recombinant Escherichia coli Outer membrane protein C protein (22-367AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein C (Outer membrane protein 1B) (Porin OmpC), ompC, meoA par
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	ompC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.