



virE2 Antibody

Product Code	CSB-PA357288LA01AYS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08062
Immunogen	Recombinant Agrobacterium fabrum Single-strand DNA-binding protein (1-556AA)
Raised In	Rabbit
Species Reactivity	Agrobacterium fabrum
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-strand DNA-binding protein (63.5 kDa virulence protein), virE2
Immunogen Species	Agrobacterium fabrum (strain C58 / ATCC 33970) (Agrobacterium tumefaciens (strain C58))
Research Area	Others
Target Names	virE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.