



# Ribosome-inactivating protein alpha-trichosanthin Antibody

<b>Product Code</b>	CSB-PA357829HA01TIF
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09989
<b>Immunogen</b>	Recombinant Trichosanthes kirilowii Ribosome-inactivating protein alpha-trichosanthin protein (24-270AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Trichosanthes kirilowii
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ribosome-inactivating protein alpha-trichosanthin (Alpha-TCS) (EC 3.2.2.22) (rRNA N-glycosidase)
<b>Immunogen Species</b>	Trichosanthes kirilowii (Chinese snake gourd) (Chinese cucumber)
<b>Research Area</b>	Others
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.