



Ribosome-inactivating protein alpha-trichosanthin Antibody, Biotin conjugated

Product Code	CSB-PA357829HD01TIF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09989
Immunogen	Recombinant Trichosanthes kirilowii Ribosome-inactivating protein alpha-trichosanthin protein (24-270AA)
Raised In	Rabbit
Species Reactivity	Trichosanthes kirilowii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein alpha-trichosanthin (Alpha-TCS) (EC 3.2.2.22) (rRNA N-glycosidase)
Immunogen Species	Trichosanthes kirilowii (Chinese snake gourd) (Chinese cucumber)
Research Area	Others