



# GPC Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA357830LC01LKW
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09991
<b>Immunogen</b>	Recombinant Lymphocytic choriomeningitis virus Pre-glycoprotein polyprotein GP complex protein (266-498AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Lymphocytic choriomeningitis virus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pre-glycoprotein polyprotein GP complex (Pre-GP-C) [Cleaved into: Stable signal peptide (SSP); Glycoprotein G1 (GP1); Glycoprotein G2 (GP2)], GPC, GP-C
<b>Immunogen Species</b>	Lymphocytic choriomeningitis virus (strain Armstrong) (LCMV)
<b>Research Area</b>	Others
<b>Target Names</b>	GPC
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.