



GPC Antibody, Biotin conjugated

Product Code	CSB-PA357830LD01LKW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09991
Immunogen	Recombinant Lymphocytic choriomeningitis virus Pre-glycoprotein polyprotein GP complex protein (266-498AA)
Raised In	Rabbit
Species Reactivity	Lymphocytic choriomeningitis virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-glycoprotein polyprotein GP complex (Pre-GP-C) [Cleaved into: Stable signal peptide (SSP); Glycoprotein G1 (GP1); Glycoprotein G2 (GP2)], GPC, GP-C
Immunogen Species	Lymphocytic choriomeningitis virus (strain Armstrong) (LCMV)
Research Area	Others
Target Names	GPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.